


new



Kemet Diamond Grinding and Fine Grinding

- Suitable for magnetic mounting

Kemet is pleased to announce the introduction of a new range of diamond grinding and fine grinding/ polishing discs for use with its metallographic grinding and polishing machines. The fine grinding and diamond polishing discs are used in combination with Kemet liquid diamond.

Benefits

- High and consistent stock removal while achieving perfect flatness.
- Maintain flatness while providing excellent surface finish.
- Excellent durability and superb value for money.
- Suitable for magnetic mounting

Diameter

100mm, 250mm & 300mm

Kemet Diamond Grinding Discs

Description	Colour	Code
Diamond Grinding Disc 200mm dia, Grit 80, 251 Micron	red	343580
Diamond Grinding Disc 250mm dia, Grit 80, 251 Micron	red	343576
Diamond Grinding Disc 300mm dia, Grit 80, 251 Micron	red	343572
Diamond Grinding Disc 200mm dia, Grit 120, 107 Micron	yellow	343581
Diamond Grinding Disc 250mm dia, Grit 120, 107 Micron	yellow	343577
Diamond Grinding Disc 300mm dia, Grit 120, 107 Micron	yellow	343573
Diamond Grinding Disc 200mm dia, Grit 220, 76 Micron	green	343582
Diamond Grinding Disc 250mm dia, Grit 220, 76 Micron	green	343578
Diamond Grinding Disc 300mm dia, Grit 220, 76 Micron	green	343574
Diamond Grinding Disc 200mm dia, Grit 600, 30 Micron	blue/black	343583
Diamond Grinding Disc 250mm dia, Grit 600, 30 Micron	blue/black	343579
Diamond Grinding Disc 300mm dia, Grit 600, 30 Micron	blue/black	343575
Diamond Grinding Disc 200mm dia, Grit 1200, 15 Micron	black	343588
Diamond Grinding Disc 250mm dia, Grit 1200, 15 Micron	black	343589
Diamond Grinding Disc 300mm dia, Grit 1200, 15 Micron	black	343590

Kemet Fine Grinding Discs

Description	Colour	Code
Fine Grinding Disc 200mm dia. Use with 15-9µm	red	343584
Fine Grinding Disc 250mm dia. Use with 15-9µm	red	343586
Fine Grinding Disc 300mm dia. Use with 15-9µm	red	343570
Fine Grinding Disc 200mm dia. Use with 6-3µm	green	343585
Fine Grinding Disc 250mm dia. Use with 6-3µm	green	343587
Fine Grinding Disc 300mm dia. Use with 6-3µm	green	343571