

## Geological Sample Holders



Geological samples holders are available for the Forcimat-TS automatic specimen mover. Each specimen holder can hold 1 slide, 27mm x 46mm.

### Accessories

- Forcimat specimen mover, 3 x 60mm diameter
- Forcimat specimen mover, 6 x 40mm diameter
- Set of insert rings for 30mm
- Set of insert rings for 25mm
- LAP-TS special specimen holder with boron carbide stops for lapping thin section. 60mm diameter
- POL-TS special specimen holder for polishing of thin sections. 60mm diameter